

EM-Tec Disc Type Sputter Target Selection Table

Manufacturer	Models	Ø60	Ø58	Ø57	Ø54	Ø50	Ø24	Ø19
Agar	Manual, Automatic			X				
Anatech/Technics	Hummer I/II/III/Jr					X		
Bal-Tec	SCD 005 / 040/050/500				X			
Bal-Tec	MED 10 / 20				X			
Bio-Rad	SC5200 / SC502			X				
Bio-Rad	E7400 / LT7400						X	
Cressington	108manual / 108auto			X				
Cressington	108auto/SE			X				
Cressington	208HR			X				
Denton	Desk II/II/IV/V	X						
Edwards	150B / Scancoat	X						
Emitech	SC7620 / SC7610			X				
Emitech on Ø60mm target holder	K500X / K550X/ K650X	X						
Emitech	K575X&XD K675X&XD			X				
Emitech	Q150			X				
EMS	150R / T 300R / T			X				
EMS	K575X&XD K675X&XD			X				
EMS	SC7620 / SC7610			X				
EMS on Ø60mm target holder	K500X / K550X / K650X	X						
Emscope on Ø58mm target holder	SC500 / SC500A		X					
Emscope	SC502			X				
Fullam	EffaCoater / EMS-76					X		
Gatan	681							X
Gatan	682 PECS							X
Gatan	685 PECS II							X
HHV	BT150 / 300				X			
HoYeon	HC21			X				
ISI	ISI 5400			X				
JEOL	JFC- 1200/1300/ 1400/1600			X				

JEOL	JFC-2300HR			X				
Leica	EM SCD 005/050/500				X			
Leica	EM MED 020				X			
Leica	ACE 200 / ACE 600				X			
Pelco	SC4 / SC5 SC6 / SC7			X				
Polaron / Bio-Rad	E5000 / E5200 / E5400 / E5800			X				
Polaron / Bio-Rad	PS 3			X				
Polaron / Quorum	SC7620 / SC7610			X				
Quorum	Q150R / T Q300R / T SC7620			X				
Quorum	PP2000 / PP2000T						X	
Quorum	PP3000T / PP3010T						X	
Safematic	CCU-010				X			
SEC	MCM-100 / MCM-200					X		
SPI	SPI-Module			X				
SPI	SPI Super					X		
VG Microtech / Fisons / Thermo	PP2200 / PP2000T						X	

Note 1: Polaron coaters have also been sold under the Bio-Rad, Fisons, VG Microtech, Thermo VG, Emscope and Quorum brand.

Note 2: Thickness of the targets can vary between 0.1 and 3.2mm